

# Nonin Medical Comparison:

## NoninConnect™ Model 3230 versus Nonin Onyx® II, Model 9560



	<b>NoninConnect™ 3230</b> with <i>Bluetooth®</i> Smart Wireless Technology	<b>Onyx® II, Model 9560</b> with <i>Bluetooth®</i> Wireless Technology
<b>Typical Application</b>	Spot check readings	Spot check readings
<b>Sensors</b>	Fixed (finger clip)	Fixed (finger clip)
<b>Wireless Radio</b>	Low energy Bluetooth (BLE 4.0)	Classic Bluetooth (2.0)
<b>Display</b>	Color Display	Monochrome Display
<b>Battery Life</b>	2200 spot check readings	600 spot check readings
<b>Warranty</b>	2-year warranty	4-year warranty
<b>Display Orientation</b>	Readings facing patient	Readings facing clinician
<b>Nonin Correct Check™</b>	Yes	No
<b>Simplified Pairing</b>	Yes	No
<b>Nonin SmartPoint™</b>	Yes	Yes
<b>Pleth Waveform Output</b>	No	Yes